



Cystatin C Calibrators KIT

✓ REFERENCE



calibrators KIT	CYREK-000	5 x 1 ml	2-8°C
Human Cystatin C in synthetic biological fluid standardized from the international standard ERM-DA471/IFCC ¹ , sodium azide (< 1g/l)			
Lot #	17H22		
Expiry date	08/2018		
Control date	03/10/2017		
Quality control report #	DGM-QAC-REP-17155		
Document prepared and signed by	Gilbert Ameryckx		

✓ SAMPLES AND REFERENCE VALUES

See the corresponding reagents technical sheet.

✓ COMPOSITION

Cystatin C calibrators are synthetic biological fluids containing human Cystatin C at fixed value diluted in HEPES pH 7.4 buffer containing stabilizers and sodium azide (<1g/l) as preservative.

✓ PRINCIPLE OF TEST

The human serum sample reacts upon colloidal gold coated with a mixture of polyclonal antibodies to Cystatin C. In the presence of cystatin C, the particles agglutinate, inducing a red shift in the absorption spectrum of the colloid. This induces an increase of the optical density at 600 nm, which is directly proportional to the ferritin concentration in the sample.

✓ PRECAUTIONS

For in vitro single diagnostic use. To be handled by entitled Personnel. Products from human source were tested and found free from HBsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective. Products containing sodium azide have to be handled with care; avoid ingestion and contact with skin and mucous membranes. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides.

✓ ANALYTICAL PERFORMANCES

See the corresponding reagents technical sheet.

✓ PREPARATION AND REAGENTS STABILITY

The calibrators are ready for use; once opened, they are stable until expiry date if stored stoppered in appropriate temperature conditions and without any contamination.

✓ METHOD OF ANALYSIS AND CALCULATION

See the corresponding reagents technical sheet.

✓ QUALITY CONTROL

Accuracy and reproducibility: Analytical performances can be checked with the internal quality control serum of the laboratory.

Calibration: Calibration curve and stability of calibration curve can be validated with the DiAgam calibration controls (CYCOS-002, CYCON-002, CYCOX-002). Calibrate when the quality control results are outside acceptable range (contact the manufacturer if the deviations persist), when the reagent lot number changes or when government regulations require.

✓ BIBLIOGRAPHY

- (1) I. Zegers et al. Certification report – Certification of cystatin C in the human serum reference material ERM®-DA471/IFCC
http://www.erm-crm.org/ERM_products/search/reports/DA471.pdf



CYREKFTEN 03/10/2017 v00

Cystatin C	ROCHE® VITROS® SYSTEMS	CAL 1	CAL 2	CAL 3	CAL 4	CAL 5
		mg/l	mg/l	mg/l	mg/l	mg/l
		0.641	1.284	1.964	2.976	5.100

Value assigned from the Cystatin C international standard ERM-DA471/IFCC.
Please contact your supplier in case of other analyzers utilization.