



CATALOGUE 2019

ALINITY LINE





Company profile

MISSION STATEMENT

«To have a perfect control of Development, Manufacturing process & Selling of In Vitro Diagnostic liquid reagents with a strong focus on our Customer vision»

Your R&D Partner

DiAgam manufactures **reagents for immunoturbidimetry**.

Our company specializes in high technology R&D, based on immunology technologies. DiAgam is the creator of liquid laboratory reagents containing **gold particles**.

We are currently taking part in a project to develop **new nanotechnologies** based on the use of **colloidal particles** in the field of medical diagnosis. This activity has also succeeded in completing our products range with Cystatin C, Ferritin & Rheumatoid factor.

Manufacturer of ready-to-use reagents

All our products are developed and manufactured within the company.

Our reagents, controls and calibrators are liquid, stable and ready to use on **all opened biochemistry ana-lyzers even with barcodes** (Abbott - Alinity – Architect , Beckman Coulter – AU series, IDS-iSYS Multi-Discipline Automated System, Menarini – Zenit RA, Mindray BS series, Ortho - VITROS® 5,1 FS Chemistry, VITROS® 4600 Chemistry, VITROS® 5600 Integrated Systems, Roche - Modular®- COBAS C, Siemens - Advia® 1800 - Advia® 2400 - Atellica®, Thermo-Konelab®, etc...)

DiAgam is a member of the National Association of **European Diagnostic Manufacturers Association (EDMA)**. All reagents are labelled and DiAgam is certified **ISO 9001:2015** and **ISO 13485:2016**.

All calibrators and controls are traceable to the **International Standards** (WHO, IFCC, NIBSC, etc.)

Worldwide reliability

DiAgam is present all around the world. We sell directly in Belgium, France, Spain & Portugal. We also export to French Speaking African countries, Argentina, Brazil, Czech Republic, Greece, Hungary, India, Mexico, Morocco, Poland, Sri Lanka, Ukraine, Vietnam, etc.

Thanks to our **manufac-turing flexibility**, we can propose all our parameters in three different packagings: **Kits** especially dedicated to particular analyzers, **Vials** filled and **Bulk**, available under **DiAgam** label or in **OEM**.

We respect our philosophy by offering our cus-tomer, a **high quality of our products** with a most **competitive price** along with a continuous quest to **satisfy their needs**.

CATALOGUE 2019

ALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ALTUR-B00/ALI	1 x 74 ml R1 + 1 x 19,5 ml R2	488	36 months	serum or plasma	1,096 to 120,44 (g/L)
	ALTUR-H00/ALI	5 x 74 ml R1 + 5 x 19,5 ml R2	2438			
	ALTUR-L00/ALI	10 x 74 ml R1 + 10 x 19,5 ml R2	4875			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ATTUR-B00/ALI	1 x 74 ml R1 + 1 x 14,5 ml R2	290	36 months	serum or plasma	0,007 to 3,5 (g/L)
	ATTUR-H00/ALI	5 x 74 ml R1 + 5 x 14,5 ml R2	1450			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	A1DGX-B00/ALI	1 x 25 ml R1 + 1 x 25 ml R1	250	6 months	serum or plasma	1,29 to 166,6 (mg/L)
	A1DGX-H00/ALI	5 x 25 ml R1 + 5 x 25 ml R2	1250			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	6 months
-----------	-----------	--------	----------

CONTROLS

M	A1CON-002	1x2 ml	6 months
L	A1COS-002	1x2 ml	6 months

CATALOGUE 2019

ALPHA-2 MACROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MGTUR-B00/ALI	1 x 74 ml R1 + 1 x 14,5 ml R2	290	36 months	serum or plasma	0,026 to 4,96 (g/L)
	MGTUR-H00/ALI	5 x 74 ml R1 + 5 x 14,5 ml R2	1450			
	MGTUR-L00/ALI	10 x 74 ml R1 + 10 x 14,5 ml R2	2900			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	12 months
URINARY HIGH	MGCU2-002	1x2 ml	12 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2			
KIT	IGLCR-B00/ALI	1 x 38 ml R1 + 1 x 8 ml R2	229	36 months	serum or plasma	3,375 to 249,78 (mg/L)
	IGLCR-H00/ALI	5 x 38 ml R1 + 5 x 8 ml R2	1143			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2			
KIT	IGLCR-B00/ALI	1 x 74 ml R1 + 1 x 13,5 ml R2	540	36 months	serum or plasma	3,375 to 249,78 (mg/L)
	IGLCR-H00/ALI	5 x 74 ml R1 + 5 x 13,5 ml R2	2700			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

CATALOGUE 2019

PREALBUMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	PATUR-B00/ALI	1 x 74 ml R1 + 1 x 14,5 ml R2		290	36 months	serum or plasma	6,53 to 1000 (mg/L)
	PATUR-H00/ALI	5 x 74 ml R1 + 5 x 14,5 ml R2		1450			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RBTUR-B00/ALI	1 x 31 ml R1 + 1 x 10 ml R2		200	24 months	serum or plasma	7,81 to 141 (mg/L)
	RBTUR-H00/ALI	5 x 31 ml R1 + 5 x 10 ml R2		1000			
	RBTUR-L00/ALI	10 x 31 ml R1 + 10 x 10 ml R2		2000			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	24 months
M	RBCON-002	1x2 ml	24 months
L	RBCOS-002	1x2 ml	24 months



Your expert in specific proteins

-
- (*) Tests :** Amount of theoretical tests. May vary from a analyzer to another one.
- (**) Shelf Life :** DiAgam provides to distributors a minimum of **12 months** .(except for Alpha1-Microglobulin). The time period written in this document shows the maximum shelf life time after manufacturing.
- (***) Sample :** DiAgam may provide some applications for other samples, such as LCR, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.
- (****) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.