



CATALOGUE 2019

ARCHITECT Line





Company profile

MISSION STATEMENT

«To have a perfect control of Development, Manufacturing process & Selling of In Vitro Diagnostic liquid reagents with a strong focus on our Customer vision»

Your R&D Partner

DiAgam manufactures **reagents for immunoturbidimetry**.

Our company specializes in high technology R&D, based on immunology technologies. DiAgam is the creator of liquid laboratory reagents containing **gold particles**.

We are currently taking part in a project to develop **new nanotechnologies** based on the use of **colloidal particles** in the field of medical diagnosis. This activity has also succeeded in completing our products range with Cystatin C, Ferritin & Rheumatoid factor.

Manufacturer of ready-to-use reagents

All our products are developed and manufactured within the company.

Our reagents, controls and calibrators are liquid, stable and ready to use on **all opened biochemistry ana-lyzers even with barcodes** (Abbott - Alinity – Architect , Beckman Coulter – AU series, IDS-iSYS Multi-Discipline Automated System, Menarini – Zenit RA, Mindray BS series, Ortho - VITROS® 5,1 FS Chemistry, VITROS® 4600 Chemistry, VITROS® 5600 Integrated Systems, Roche - Modular®-COBAS C, Siemens - Advia® 1800 - Advia® 2400 - Atellica®, Thermo-Konelab®, etc...)

DiAgam is a member of the National Association of **European Diagnostic Manufacturers Association (EDMA)**. All reagents are labelled and DiAgam is certified **ISO 9001:2015** and **ISO 13485:2016**.

All calibrators and controls are traceable to the **International Standards** (WHO, IFCC, NIBSC, etc.)

Worldwide reliability

DiAgam is present all around the world. We sell directly in Belgium, France, Spain & Portugal. We also export to French Speaking African countries, Argentina, Brazil, Czech Republic, Greece, Hungary, India, Mexico, Morocco, Poland, Sri Lanka, Ukraine, Vietnam, etc.

Thanks to our **manufac-turing flexibility**, we can propose all our parameters in three different packagings: **Kits** especially dedicated to particular analyzers, **Vials** filled and **Bulk**, available under **DiAgam** label or in **OEM**.

We respect our philosophy by offering our cus-tomer, a **high quality of our products** with a most **competitive price** along with a continuous quest to **satisfy their needs**.

ALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	ALTUR-B00/ARC	1 x 40 ml R1 + 1 x 10 ml R2	250	36 months	serum or plasma	0.40 to 121.71 (g/l)
	ALTUR-H00/ARC	1 x 80 ml R1 + 1 x 20 ml R2	500			
	ALTUR-L00/ARC	7 x 80 ml R1 + 7 x 20 ml R2	3500			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOÏD

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	GPTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2	150	36 months	serum	0.143 to 2.42 (g/l)
	GPTUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2	600			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	ATTUR-B00/ARC	1 x 45 ml R1 + 1 x 9 ml R2	180	24 months	serum or plasma	0.007 to 3.5 (g/l)
	ATTUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	A1DGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2	150	6 months	serum or plasma	1,29 to 166,6 (mg/l)
	A1DGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2	550			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	6 months
-----------	-----------	--------	----------

CONTROLS

M	A1CON-002	1x2 ml	6 months
L	A1COS-002	1x2 ml	6 months

ALPHA-2 MACROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MGTUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2	220	36 months	serum or plasma	0.026 to 4.96 (g/l)
	MGTUR-L00/ARC	10 x 55 ml R1 + 10 x 11 ml R2	2200			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	12 months
URINARY HIGH	MGCU2-002	1x2 ml	12 months

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ASLTX-B00/ARC	1 x 70 ml R1 + 1 x 8 ml R2	296	20 months	serum or plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/ARC	2 x 87 ml R1 + 2 x 10 ml R2	741			

CALIBRATORS

H	ASREH-001(-005)	1x1 or (1x5) ml	12 months
----------	-----------------	-----------------	-----------

CONTROLS

M	ASCON-002	1x2 ml	12 months
----------	-----------	--------	-----------

ANTITHROMBIN III

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2			
KIT	A3TUR-B00/ARC	1 x 50 ml R1 + 1 x 10 ml R2	250	36 months	plasma	28.15 to 806.30 (mg/l)
	A3TUR-H00/ARC	2 x 70 ml R1 + 2 x 15 ml R2	700			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	AATUR-B00/ARC	1 x 38 ml R1 + 1 x 10 ml R2	152	36 months	serum or plasma	0.029 to 400 (mg/dl)
	AATUR-H00/ARC	2 x 75 ml R1 + 2 x 20 ml R2	600			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ABTUR-B00/ARC	1 x 41 ml R1 + 1 x 8 ml R2	242	36 months	serum or plasma	0.019 to 4.0 (g/l)
	ABTUR-H00/ARC	2 x 67 ml R1 + 2 x 13 ml R2	788			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

BETA 2 MICROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	B2DGX-B00/ARC	1 x 27 ml R1 + 1 x 12 ml R2	150	15 months	serum or plasma	0,132 to 95 (mg/l)
	B2DGX-H00/ARC	2 x 45 ml R1 + 2 x 20 ml R2	500			
	B2DGX-L00/ARC	10 x 45 ml R1 + 10 x 20 ml R2	2500			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	18 months
M	B2CON-003	1x3 ml	18 months
L	B2COS-003	1x3 ml	18 months

CATALOGUE 2019

BETA 2 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	BUDGX-B00/ARC	R1 1 x 23 ml R1 + 1 x 10 ml R2	167	15 months	Urine	0.055 to 9.3 (mg/l)
	BUDGX-H00/ARC	R2 2 x 23 ml R1 + 2 x 10 ml R2	333			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	18 months
L	B2COS-003	1x3 ml	18 months

CERULOPLASMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	CETUR-B00/ARC	R1 1 x 40 ml R1 + 1 x 8 ml R2	200	36 months	serum	9.48 to 1051.20 (mg/l)
	CETUR-H00/ARC	R2 2 x 60 ml R1 + 2 x 12 ml R2	600			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	C3TUR-B00/ARC	R1 1 x 41 ml R1 + 1 x 8 ml R2	216	36 months	serum or plasma	0.015 to 4.31 (g/l)
	C3TUR-H00/ARC	R2 2 x 56 ml R1 + 2 x 11 ml R2	595			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample**	Analytical Range
KIT		<i>R1</i>	<i>R2</i>			
	C4TUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2	200	36 months	serum or plasma	0.002 to 0.825 (g/l)
	C4TUR-H00/ARC	2 x 60 ml R1 + 2 x 12 ml R2	600			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CPDGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2	150	16 months	serum or plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2	550			
CPDGX-L00/ARC	10 x 55 ml R1 + 10 x 55 ml R2	5500				

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	18 months
M	CPCON-002	1x2 ml	18 months
L	CXCON-002	1x2 ml	18 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CPTXL-B00/ARC	1 x 55,5 ml R1 + 1 x 8 ml R2	222	36 months	serum or plasma	0.607 to 425 (mg/l)
	CPTXL-H00/ARC	2 x 85 ml R1 + 2 x 12 ml R2	667			
CPTXL-L00/ARC	7 x 85 ml R1 + 7 x 12 ml R2	2333				

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	18 months
M	CPCON-002	1x2 ml	18 months
L	CXCON-002	1x2 ml	18 months

CATALOGUE 2019

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i> <i>R2</i>				
	CYCOL-B00/ARC	1 x 27 ml R1 + 1 x 11 ml R2	147	6 months	serum or plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/ARC	2 x 48 ml R1 + 2 x 20 ml R2	533			
CYCOL-L00/ARC	10 x 48 ml R1 + 10 x 20 ml R2	2667				

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	12 months
------------------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	12 months
M	CYCON-002	1x2 ml	12 months
L	CYCOS-002	1x2 ml	12 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample**	Analytical Range
KIT		<i>R1</i> <i>R2</i>				
	FNCOL-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2	200	12 months	serum	6.344 to 500 (µg/l)
FNCOL-H00/ARC	2 x 55 ml R1 + 2 x 11 ml R2	550				

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
----------	-----------	--------	-----------

FIBRINOGEN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i> <i>R2</i>				
	FITUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2	150	36 months	plasma (from sodium citrate)	0.05 to 10 (g/l)
FITUR-H00/ARC	2 x 50 ml R1 + 1 x 15 ml R2	300				

CALIBRATORS

M	FIREM-001	1x1 ml	12 months
----------	-----------	--------	-----------

CONTROLS

M	FICON-001	1x1 ml	12 months
----------	-----------	--------	-----------

CATALOGUE 2019

HAPTOGLOBIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	HATUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2		150	36 months	serum	0.002 to 3.86 (g/l)
	HATUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2		600			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample**	Analytical Range
		R1	R2				
KIT	IATUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2		229	36 months	serum	0.02 to 8.0 (g/l)
	IATUR-H00/ARC	2 x 54 ml R1 + 2 x 11 ml R2		629			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample**	Analytical Range
		R1	R2				
KIT	IGTUR-B00/ARC	1 x 25 ml R1 + 1 x 9,5 ml R2		148	36 months	serum	0.125 to 35 (g/l)
	IGTUR-H00/ARC	2 x 50 ml R1 + 2 x 19 ml R2		594			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2			
KIT	IGLCR-B00/ARC	1 x 74 ml R1 + 1 x 13,5 ml R2	540	36 months	serum	3,375 to 249,78 (mg/l)
	IGLCR-H00/ARC	5 x 74 ml R1 + 5 x 13,5 ml R2	2700			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2			
KIT	IMTUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2	229	36 months	serum	0.038 to 7.0 (g/l)
	IMTUR-H00/ARC	2 x 49,5 ml R1 + 2 x 10,5 ml R2	600			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	LPTUR2-B00/ARC	1 x 64 ml R1 + 1 x 8 ml R2	320	12 months	serum or plasma	2.35 to 94.60 (mg/dl)
	LPTUR2-H00/ARC	2 x 80 ml R1 + 2 x 10 ml R2	800			

CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	10 months
------------------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	12 months
----------	-----------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MATUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2	229	36 months	urine	1.86 to 400 (mg/l)
	MATUR-H00/ARC	2 x 52,5 ml R1 + 2 x 10,5 ml R2	600			
	MATUR-L00/ARC	7 x 75 ml R1 + 7 x 15 ml R2	3000			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	18 months
M	MACON-002	1x2 ml	18 months
L	MACOS-002	1x2 ml	18 months

CATALOGUE 2019

PREALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	PATUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2	220	36 months	serum and plasma	6.53 to 1000 (mg/l)

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml		24 months		
-----------	-----------	--------	--	-----------	--	--

CONTROLS

H	MPCOX-002	1x2 ml		24 months		
M	MPCON-002	1x2 ml		24 months		
L	MPCOS-002	1x2 ml		24 months		

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	RBTUR-B00/ARC	1 x 31 ml R1 + 1 x 10 ml R2	200	24 months	urine	7,81 to 141 (mg/l)
	RBTUR-H00/ARC	2 x 31 ml R1 + 2 x 10 ml R2	400			
	RBTUR-L00/ARC	10 x 31 ml R1 + 10 x 10 ml R2	2000			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml		24 months		
-----------	-----------	--------	--	-----------	--	--

CONTROLS

H	RBCOX-002	1x2 ml		24 months		
M	RBCON-002	1x2 ml		24 months		
L	RBCOS-002	1x2 ml		24 months		

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	RFCOL-B00/ARC	1 x 30 ml R1 + 1 x 20 ml R2	150	18 months	serum and plasma	2.228 to 120 (kUI/l)
	RFCOL-H00/ARC	2 x 60 ml R1 + 2 x 40 ml R2	600			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml		24 months		
-----------	-----------	--------	--	-----------	--	--

CONTROLS

M	RFCON-002	1x2 ml		24 months		
---	-----------	--------	--	-----------	--	--

TRANSFERRIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample**	Analytical Range
		<i>R1</i>	<i>R2</i>			
KIT	TFTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 10,5 ml R2	150	36 months	serum and plasma	0.076 to 8.0 (g/l)
	TFTUR-H00/ARC	2 x 72 ml R1 + 2 x 20 ml R2	571			
	TFTUR-L00/ARC	7 x 72 ml R1 + 7 x 20 ml R2	2000			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml		24 months		
-----------	-----------	--------	--	-----------	--	--

CONTROLS

H	MPCOX-002	1x2 ml		24 months		
M	MPCON-002	1x2 ml		24 months		
L	MPCOS-002	1x2 ml		24 months		



Your expert in specific proteins

(*) Tests : Amount of theoretical tests. May vary from a analyzer to another one.

() Shelf Life :** DiAgam provides to distributors a minimum of **12 months** (except for Cystatin-c, Ferritin and Lp(a): < 12 months). The time period written in this document shows the maximum shelf life time after manufacturing.

(*) Sample :** DiAgam may provide some applications for other samples, such as LCR, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.

(**)**
Analytical range : For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.

IVD Medical device for in vitro diagnosis.



DIAGAM SA, Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.