



*Your expert in specific proteins*

## CATALOGUE 2019



*Reagents, Calibrators & Controls*

Immunoturbidimetry

**BULK OEM**



## ALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	17BUF-999	Bulk	2-25 °C	36 months	serum / plasma	1.096 to 119.79 (g/l)	150/40 (µl)
<b>Ab</b>	ALTUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## ALPHA1-ANTITRYPSIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.007 to 3.5 (g/l)	250/50 (µl)
<b>Ab</b>	ATTUR-H99	Bulk	2-8 °C	24 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## ALPHA1-GLYCOPROTEIN ACID / OROSOMUCOID

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	GYBUF-999	Bulk	2-25 °C	36 months	serum	0.143 to 2.0 (g/l)	250/50 (µl)
<b>Ab</b>	GPTUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------



### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## ALPHA-1 MICROGLOBULIN (URINE)

**Contact us**

## ALPHA2-MACROGLOBULIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.026 to 4.96 (g/l)	250/50 (µl)
<b>Ab</b>	MGTUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months
URINARY LOW	MGCU1-999	Bulk	2-8 °C	12 months
URINARY HIGH	MGCU2-999	Bulk	2-8 °C	12 months

## ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	ASBUF-999	Bulk	2-8 °C	20 months	serum	10 to 900 (IU/ml)	235/30 (µl)
<b>Ab</b>	ASLTX-H99	Bulk	2-8 °C	20 months			



### CALIBRATORS

H	ASREH-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

### CONTROLS

M	ASCON-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

## ANTITHROMBIN III

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TIBUF-999	Bulk	2-25 °C	36 months	plasma	28.15 to 806.3 (mg/l)	250/50 (µl)
<b>Ab</b>	A3TUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

# CATALOGUE 2019 - Bulk OEM

## APOLIPOPROTEIN A1

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.029 to 400 (mg/dl)	250/66 (µl)
<b>Ab</b>	AATUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

## APOLIPOPROTEIN B

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.019 to 4.0 (g/l)	250/50 (µl)
<b>Ab</b>	ABTUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

## BETA2-MICROGLOBULIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	B2BUX-999	Bulk	2-8 °C	15 months	serum / plasma	0,132 to 95 (mg/l)	180/80 (µl)
<b>Ab</b>	B2DGX- H99	Bulk	2-8 °C	15 months			



### CALIBRATORS

KIT 5 LEV	B2REK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------



### CONTROLS

H	B2COX-001	lyophilized	2-8 °C	18 months
M	B2CON-003	lyophilized	2-8 °C	18 months
L	B2COS-003	lyophilized	2-8 °C	18 months

## BETA2-MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	BUBUX-999				Urine	0,055 - 9,3 mg/L	135/60 (µl)
<b>Ab</b>	BUDGX- H99			Contact us			

## CERULOPLASMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	MEBUF-999	Bulk	2-8 °C	36 months	serum	9.48 to 1051.2 (mg/l)	250/50 (µl)
<b>Ab</b>	CETUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## COMPLEMENT C3

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7BUF-999	Bulk	2-25 °C	36 months	serum / plasma	0.015 to 4.31 (g/l)	250/50 (µl)
<b>Ab</b>	C3TUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## COMPLEMENT C4

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7BUF-999	Bulk	2-25 °C	36 months	serum / plasma	0.002 to 0.825 (g/l)	250/50 (µl)
<b>Ab</b>	C4TUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

# CATALOGUE 2019 - Bulk OEM

## CRP GOLD LATEX

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	CPBUF-999	Bulk	2-8 °C	16 months	serum / plasma	0.055 to 290 (mg/l)	100/100 (µl)
<b>Ab</b>	CPDGX-H99	Bulk	2-8 °C	16 months			

### CALIBRATORS

KIT 5 LEV	CPRLK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCXN-999	Bulk	2-8 °C	18 months

## CRP XL WIDE RANGE

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	PITRX-999	Bulk	2-25 °C	36 months	serum / plasma	0.607 to 300 (mg/l)	250/36 (µl)
<b>Ab</b>	CPTXL-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	CPREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCXN-999	Bulk	2-8 °C	18 months

## CRP XS ULTRA - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	PITRX-999	Bulk	2-25 °C	36 months	serum / plasma	0.471 to 285.05 (mg/l)	250/20 (µl)
<b>Ab</b>	CPTXS-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS



KIT 5 LEV	CPR(EK/LK)-999	Bulk	2-8 °C	18 months
-----------	----------------	------	--------	-----------

### CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCXN-999	Bulk	2-8 °C	18 months

# CATALOGUE 2019 - Bulk OEM

## CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum	0.15 to 9 (mg/l)	180/75 (µl)
<b>Ab</b>	CYCOL-H99	Bulk	2-8 °C	6 months			



### CALIBRATORS

KIT 5 LEV	CYREK-999	Bulk	2-8 °C	12 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	CYCOX-999	Bulk	2-8 °C	12 months
M	CYCON-999	Bulk	2-8 °C	12 months
L	CYCOS-999	Bulk	2-8 °C	12 months

## FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	FDBUF-999	Bulk	2-8 °C	18 months	serum	6.614 to 500 (µg/l)	250/50 (µl)
<b>Ab</b>	FNCOL-H99	Bulk	2-8 °C	12 months			



### CALIBRATORS

KIT 5 LEV	FNRGK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

M	FNCON-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

## FIBRINOGEN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	plasma	0.5 to 10 (g/l)	250/80 (µl)
<b>Ab</b>	FITUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

M	FIREM-001	lyophilized	2-8 °C	12 months
---	-----------	-------------	--------	-----------

### CONTROLS

M	FICON-001	lyophilized	2-8 °C	12 months
---	-----------	-------------	--------	-----------

## HAPTOGLOBIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7TRX-999	Bulk	2-25 °C	36 months	serum	0.002 to 4.0 (g/l)	250/50 (µl)
<b>Ab</b>	HATUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## IMMUNOGLOBULINS A

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.02 to 8.0 (g/l)	250/50 (µl)
<b>Ab</b>	IATUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## IMMUNOGLOBULINS G

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.125 to 35 (g/l)	250/95 (µl)
<b>Ab</b>	IGTUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------



### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

## IMMUNOGLOBULINS G (CSF)

**Contact us**

## IMMUNOGLOBULINS M

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.038 to 7.0 (g/l)	250/50 (µl)
<b>Ab</b>	IMTUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------



### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months



## CATALOGUE 2019 - Bulk OEM

### LIPOPROTEIN (A)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	2.35 to 94.60 (mg/dl)	200/25 (µl)
<b>Ab</b>	LPTUR2-H99	Bulk	2-8 °C	12 months			



### CALIBRATORS

KIT 5 LEV	LPREK-999	Bulk	2-8 °C	10 months
-----------	-----------	------	--------	-----------

### CONTROLS

M	LPCON-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

### MICROALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7BUF-999	Bulk	2-25 °C	36 months	urine	1.86 to 400 (mg/l)	250/50 (µl)
<b>Ab</b>	MATUR-H99	Bulk	2-8 °C	36 months			



### CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MACOX-999	Bulk	2-8 °C	18 months
M	MACON-999	Bulk	2-8 °C	18 months
L	MACOS-999	Bulk	2-8 °C	18 months

### PREALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7BUF-999	Bulk	2-25 °C	36 months	serum	6.53 to 1000 (mg/l)	250/50 (µl)
<b>Ab</b>	PATUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------



### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

### RETINOL BINDING PROTEIN

**Contact us**

## RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	TRBUF-999	Bulk	2-8 °C	36 months	serum	2.092 to 120 (kUI/l)	150/100 (µl)
<b>Ab</b>	RFCOL-H99	Bulk	2-8 °C	18 months			



### CALIBRATORS

KIT 5 LEV	RFREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

M	RFCON-999	Bulk	2-8 °C	24 months
---	-----------	------	--------	-----------

## TRANSFERRIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
<b>BUF</b>	I7BUF-999	Bulk	2-25 °C	36 months	serum	0.076 to 8.0 (g/l)	250/70 (µl)
<b>Ab</b>	TFTUR-H99	Bulk	2-8 °C	36 months			

### CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

### CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months



*Your expert in specific proteins*

**(\*) Shelf Life :** DiAgam provides to distributors a minimum of *12 months (except for Cystatin-c, Ferritin and Lp(a): < 12 months)* . The time period written in this document shows the maximum shelf life time after manufacturing.

**(\*\*) Sample :** DiAgam may provide some applications for other samples, such as LCR, Urin, Veterany samples etc.... Please contact your distributor or us directly at *mail@diagam.com* or dial *+32 68 55 14 82*.

**(\*\*\*) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at *mail@diagam.com* or dial *+32 68 55 14 82* in order to get the updated data.

**(\*\*\*\*)** Bulk formats or any production require a minimum volume of order of 250 mL for calibrators and controls. If less than 250 mL of calibrators and controls are ordered, DiAgam may refuse to process the order.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel, 7822 Ghislenghien BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© DiAgam. 2019. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.